

How to use these sheets

One sheet per unit operation, ten in all, covering the 21 canonical CP-2 steps (buffer prep, D0.1-D0.5, Days 1-4, D4.1-D4.15). The left half restates the step as RTx described it (clock, grade, equipment, in-process tests, the change vs CP-1). The right half is ruled space for what the room decides. Write in pen; the sheets are the workshop record and are transcribed into the action log within 24 h.

Clocks: 'RTx stated' figures come from the RTx 'US CDMO - Pivotal Process Detailed' docx (3-4 Sep 2026) and briefing pack slides 4-9. 'MADE working estimate' figures come from the MADE capacity model and are to be confirmed at the Day-4 rehearsal; they are not facts. Day 4 elapsed time: the RTx PFD implies ~20 h harvest-to-CRF; the MADE working estimate is 12 h (target 11.5-14 h).

Index

0	Buffer and reagent preparation	Day -15 to -1	RTx 3 h	MADE n/a	none
1	Starting-material assessment and CD14+ enrichment	Day 0	RTx 10.83 h	MADE n/a	none
2	Monocyte seeding into MACS1000 bags via Prodigy formulation	Day 0	RTx 1 h (TBC)	MADE n/a	C1, C2, C3
3	Monocyte-to-macrophage differentiation	Days 1-4	RTx 96 h	MADE n/a	C4
4	Harvest, pooling and Rotea wash into electroporation buffer	Day 4	RTx 7.33 h	MADE 3.33 h	none
5	Transfection preparation and Flowfect electroporation	Day 4	RTx 1.67 h	MADE 2.25 h	C5, C9
6	Recovery culture	Day 4	RTx 5.17 h (TBC)	MADE 5.17 h	C6
7	Harvest, Rotea wash into saline and dose setting	Day 4 into 5	RTx 2.83 h	MADE 3.75 h	none
8	DS and FDP formulation and fill on Ares X20	Day 4 into 5	RTx 3.17 h (TBC)	MADE 1.92 h	C1, C7
9	Visual inspection, cassette, controlled-rate freeze and VPLN2	Day 4 into 5	RTx 0.58 h	MADE 1.75 h	C8

The eight changes (C1-C8) and the mRNA supply variable (C9)

C1 · Liquid handling: Fresenius Cue for fluid transfer, formulation and fill to Cue replaced (Prodigy formulation bar at D0; Ares X20 at fill). RTx rationale: More appropriate options available; replacements being tested; data pending.

C2 · Culture bag fill: CellConnect / manual fill to CliniMACS Prodigy formulation bar. RTx rationale: Closed fill from Prodigy.

C3 · Culture bag size: up to 18-20 MACS500 to up to 10 MACS1000. RTx rationale: Reduce bag count and handling.

C4 · Macrophage conversion period: 5 days to 4 days. RTx rationale: Plateau in differentiation markers; reduced facility time.

C5 · Electroporation: Lonza 4D-Nucleofector LV to Kytopen Flowfect (TBC). RTx rationale: Faster throughput, qualification, customization.

C6 · Post-EP recovery / time to cryo: Overnight on magnetic stirrer (D5-D6) to ~4 h (2-6 h window) TBC. RTx rationale: Optimizing mRNA quantity in DP; IL-10 OOS risk at lower limit; capitalize on expression before mRNA decays.

C7 · Formulation / fill automation: Cue or manual to Xiogenix Ares X20. RTx rationale: Reduce manual handling, increase sterility assurance.

C8 · Cryopreservation: Planer CRF, current cassette to TF CryoMed LN2 CRF (or other) + metal cassette. RTx rationale: CDMO availability and broader access for pivotal.

C9 · mRNA supply (not a process change but a comparability variable): MS-1 (EU provider) to MS-2 (new provider). RTx rationale: Commercial-appropriate supply.

Sources: ruby_process.json v1.0 (MADE, 26 Sep 2026); RTx 'US CDMO - Pivotal Process Detailed' docx (S. Mellon / S. Howe, 3-4 Sep 2026); RTx Workshop Briefing Pack vFinal slides 4-9 (4 Sep 2026); MADE capacity model working estimates (Sep 2026); RTx x MADE pre-briefing (21 Sep 2026). US English. Equipment named is the CP-2 process equipment as stated by RTx.

0 Buffer and reagent preparation

A (ISO 5)

Day -15 to -1

Clock

RTx stated: 3 h**Steps (RTx docx)**

PREP · Process buffer, DS formulation buffer, electroporation buffer, base culture medium, recovery medium · 3 h

Starting material and inputs

- Raw materials (buffers, media, reagents)
- RTX001_01 mRNA vial thawed and transferred to a weldable storage pouch (Grade A; new in CP-2)

In-process test / sample

no IPC at this unit operation

Key equipment and consumable

Peristaltic pump; Balance; Tube welder; Heat sealer; Class A MSC

The change vs CP-1

no numbered change at this unit operation

Note

Buffers and media are prepared ahead of the run in a Grade A environment. New in CP-2: the mRNA vial is thawed and transferred to a weldable storage pouch so it can be welded into the closed Day 4 transfection set-up.

RTx change and data status**MADE view****Optimization ideas****Operational impact**

facility / equipment / staffing / QC

*facility**equipment**staffing**QC***Data needed****Decision****Owner****Date**

1 Starting-material assessment and CD14+ enrichment

C (ISO 7)

Day 0

Clock

RTx stated: 10.83 h

Steps (RTx docx)

- D0.1 · Equipment set-up · 0.5 h
- D0.2 · Cellular starting-material assessment · 0.5 h
- D0.3 · Pre-enrichment prep (loading, priming) · 3.33 h
- D0.4 · CD14+ enrichment run (selection + wash) · 6.5 h

Starting material and inputs

- Leukapheresis cell material: receipt within 48 h of collection; apheresis Tue/Wed (pivotal), arrival next afternoon; Day 0 start within 24 h of receipt

In-process test / sample

- CD45+ count (NC-202)
- %CD14+ (BD FACSLytic)
- CD45+ count post-enrichment

Key equipment and consumable

CliniMACS Prodigy (Miltenyi); NC-202 (ChemoMetec); BD FACSLytic; Class A MSC

The change vs CP-1

no numbered change at this unit operation

Note

Prodigy CD14 input envelope (MATCH precedent): up to 20e9 WBC, under 400e6 WBC/mL, up to 3.5e9 CD14+, 50-300 mL. Starting-material assessment sets the CD45+ count and %CD14+ before enrichment.

RTx change and data status

MADE view

Optimization ideas

Operational impact

facility / equipment / staffing / QC

facility

equipment

staffing

QC

Data needed

Decision

Owner

Date

2 Monocyte seeding into MACS1000 bags via Prodigy formulation bar

C (ISO 7)

Day 0

Clock

RTx stated: 1 h (TBC)

Steps (RTx docx)

D0.5 · Monocyte seeding into up to 10 MACS1000 differentiation bags via Prodigy formulation bar · 1 h TBC

Starting material and inputs

· Up to 10 x MACS1000 differentiation bags (replaces up to 18-20 x MACS500)

In-process test / sample

· Cell count end of Day 0 (NC-202)

Key equipment and consumable

Prodigy formulation bar (custom application from Miltenyi);
MACS1000 bags; CO2 incubator

The change vs CP-1

C1 Liquid handling: Fresenius Cue for fluid transfer, formulation and fill to Cue replaced (Prodigy formulation bar at D0; Ares X20 at fill)

RTx rationale: More appropriate options available; replacements being tested; data pending. RTx asks: CDMO experience of equipment?

C2 Culture bag fill: CellConnect / manual fill to CliniMACS Prodigy formulation bar

RTx rationale: Closed fill from Prodigy. RTx strategy: Updated CAP required, timelines, limited bag number, filtered air required. RTx asks: CDMO input: space / Prodigy usage

C3 Culture bag size: up to 18-20 MACS500 to up to 10 MACS1000

RTx rationale: Reduce bag count and handling. RTx strategy: Checking CQA impact (no issues so far). RTx asks: Incubator / manual handling preference / facility fit

Note

Replaces up to 18-20 MACS500 bags filled via Cue/CellConnect. A custom application (CAP) is required from Miltenyi for the formulation bar; filtered air is required. Bags go to the CO2 incubator.

RTx change and data status

MADE view

Optimization ideas

Operational impact

facility / equipment / staffing / QC

facility

equipment

staffing

QC

Data needed

Decision

Owner

Date

3 Monocyte-to-macrophage differentiation

C (ISO 7)

Days 1-4

Clock

RTx stated: 96 h

Steps (RTx docx)

D1-4 · 4-day culture in serum-free medium + M-CSF (was 5 days in CP-1) · 96 h

Starting material and inputs

· Serum-free medium + M-CSF (MATCH precedent: TexMACS + 100 ng/mL M-CSF, 50% feeds d2/d4). Feed schedule for the 4-day culture is not stated in the CP-2 docx: to confirm

In-process test / sample

no IPC at this unit operation

Key equipment and consumable

CO2 incubator (bags per shelf and incubators per batch: open);
MACS1000 differentiation bags

The change vs CP-1

C4 Macrophage conversion period: 5 days to 4 days

RTx rationale: Plateau in differentiation markers; reduced facility time. RTx strategy: Ongoing CQA / molecular pathway (omics) data. RTx asks: Throughput and starts

Note

4 days (was 5). RTx rationale: plateau in differentiation markers and reduced facility time; CQA and omics data ongoing. Incubator fit (bags per shelf, incubators per batch) is open and drives the facility model.

RTx change and data status

MADE view

Optimization ideas

Operational impact

facility / equipment / staffing / QC

facility

equipment

staffing

QC

Data needed

Decision

Owner

Date

4 Harvest, pooling and Rotea wash into electroporation buffer

C (ISO 7)

Day 4

Clock

RTx stated: 7.33 h*MADE working estimate: 3.33 h (rehearsal to confirm)*

Steps (RTx docx)

D4.1 · Harvest set-up · 0.33 h

D4.2 · Macrophage harvest and pooling (up to 10 bags) · 2.75 h

MADE 1.75 h (working estimate)

D4.3 · Rotea wash and buffer exchange into electroporation buffer · 4.25 h

MADE 1.25 h (working estimate)

Starting material and inputs

· Wash buffer; electroporation buffer

In-process test / sample

· Cell count start of Day 4 (NC-202)

· Cell count after buffer exchange

Key equipment and consumable

BSC / bench; Tube welder; Rotea (Thermo Fisher); NC-202

The change vs CP-1

no numbered change at this unit operation

Note

RTx comment on D4.2: the 7 h harvest-plus-wash figure is MACS500 time and is expected to be less with MACS1000. RTx comment on D4.3: the wide range likely reflects variable cell loads.

RTx change and data status

MADE view

Optimization ideas

Operational impact

facility / equipment / staffing / QC

*facility**equipment**staffing**QC*

Data needed

Decision

Owner

Date

5 Transfection preparation and Flowfect electroporation

C (ISO 7)

Day 4

Clock

RTx stated: 1.67 h*MADE working estimate: 2.25 h (rehearsal to confirm)*

Steps (RTx docx)

D4.4 · Transfection set-up and cell/mRNA/buffer preparation · 1.5 h

D4.5 · Electroporation (Kytopen Flowfect, continuous flow) · 0.17 h

MADE 0.75 h (working estimate)

Starting material and inputs

- mRNA RTX001_01 (RTx-supplied). MS-1 to MS-2 supply change is a comparability variable (C9)
- Electroporation buffer

In-process test / sample

- Cell count before transfection
- Post-transfection cell count

Key equipment and consumable

Balance; Tube welder; NC-202; Kytopen Flowfect electroporator (continuous flow)

The change vs CP-1

C5 Electroporation: Lonza 4D-Nucleofector LV to Kytopen Flowfect (TBC)

RTx rationale: Faster throughput, qualification, customization. RTx strategy: Equipment being tested; CDMO quality interactions. RTx asks: Any concerns?

C9 mRNA supply (not a process change but a comparability variable): MS-1 (EU provider) to MS-2 (new provider)

RTx rationale: Commercial-appropriate supply. RTx strategy: Small-scale comparison underway.

Note

Adjust to target concentration (patent: 50-100e6/mL small scale, 200e6/mL large-scale fluidics), load the cartridge and calculate the mRNA volume for the mRNA:cell ratio. Flowfect run ~5 min plus ~5 min close-out; RTx total transfection block 1.75 h. Replaces the Lonza 4D-Nucleofector LV.

RTx change and data status

MADE view

Optimization ideas

Operational impact

facility / equipment / staffing / QC

facility

equipment

staffing

QC

Data needed

Decision

Owner

Date

6 Recovery culture

C (ISO 7)

Day 4

Clock

RTx stated: 5.17 h (TBC)*MADE working estimate: 5.17 h (rehearsal to confirm)*

Steps (RTx docx)

D4.6 · Recovery medium addition and transfer to up to 4 MACS bags · 1.17 h

D4.7 · Recovery culture · 4 h TBC

MADE 4 h (working estimate)

Starting material and inputs

· Recovery medium; up to 4 MACS bags (how many and which size: open)

In-process test / sample

no IPC at this unit operation

Key equipment and consumable

Tube welder; Balance; CO2 incubator

The change vs CP-1

C6 Post-EP recovery / time to cryo: Overnight on magnetic stirrer (D5-D6) to ~4 h (2-6 h window) TBC

RTx rationale: Optimizing mRNA quantity in DP; IL-10 OOS risk at lower limit; capitalize on expression before mRNA decays. RTx strategy: CQA data collection. RTx asks: Batch throughput and staffing; does it fit CDMO shifts?

Note

RTx rationale for the ~4 h recovery (2-6 h window, TBC): mRNA residency decays (~3% of copies at 16 h) while IL-10 secretion is cumulative to ~16 h; use expression before decay and reduce OOS risk at the IL-10 lower limit. Bags returned to the incubator (how many, which size) is open.

RTx change and data status**MADE view****Optimization ideas****Operational impact**

facility / equipment / staffing / QC

*facility**equipment**staffing**QC***Data needed****Decision****Owner****Date**

7 Harvest, Rotea wash into saline and dose setting

C (ISO 7/8)

Day 4 into 5

Clock

RTx stated: 2.83 h*MADE working estimate: 3.75 h (rehearsal to confirm)*

Steps (RTx docx)

D4.8 · Harvest set-up · 0.5 h

D4.9 · Harvest and pooling into harvest bag · 1 h

D4.10 · Rotea wash and concentrate into saline (intermediate bag) · 0.33 h

MADE 1.25 h (working estimate)

D4.11 · Post-Rotea to loading on Ares X20 (waiting for stained TruCount) · 1 h

Starting material and inputs

· Saline (intermediate bag)

In-process test / sample

· Post-harvest count

· Post-Rotea count

· TruCount viable macrophage enumeration (BD FACSLytic): gates the fill

Key equipment and consumable

BSC / bench; Rotea (Thermo Fisher); NC-202; BD FACSLytic

The change vs CP-1

no numbered change at this unit operation

Note

Dose setting is performed on drug substance: the stained TruCount result gates the fill. D4.11 is the wait for TruCount between Rotea and loading on the Ares X20.

RTx change and data status

MADE view

Optimization ideas

Operational impact

facility / equipment / staffing / QC

*facility**equipment**staffing**QC*

Data needed

Decision

Owner

Date

8 DS and FDP formulation and fill on Ares X20

C (ISO 7/8)

Day 4 into 5

Clock

RTx stated: 3.17 h (TBC)*MADE working estimate: 1.92 h (rehearsal to confirm)*

Steps (RTx docx)

D4.12 · DS formulation (dilute in DS formulation buffer) then FDP formulation (CryoStor CS10 addition) · 2.5 h TBC

MADE 1.25 h (working estimate)

D4.13 · Container prep and Ares set-up; bag fill (bags 1-20); QC vial fill · 0.67 h

Starting material and inputs

- DS formulation buffer; CryoStor CS10
- CryoMACS 50 bags (10 mL fill at 1e7 cells/mL; 4-9 bags per RTx slide 6, up to 20 fill positions); ~20 x 1 mL QC vials

In-process test / sample

- QC samples for release and characterization taken here (see panel)
- Cell count results

Key equipment and consumable

Ares X20 (Xiogenix); Heat sealer; Balance

The change vs CP-1

- C1 Liquid handling: Fresenius Cue for fluid transfer, formulation and fill to Cue replaced (Prodigy formulation bar at D0; Ares X20 at fill)**

RTx rationale: More appropriate options available; replacements being tested; data pending. RTx asks: CDMO experience of equipment?

- C7 Formulation / fill automation: Cue or manual to Xiogenix Ares X20**

RTx rationale: Reduce manual handling, increase sterility assurance. RTx strategy: Water/buffer runs complete, cell tests underway. RTx asks: CDMO view on equipment and space

Note

DP: CryoMACS 50 bags, 10 mL fill at 1e7 cells/mL, multiple bags (4-9 per RTx slide 6, up to 20 fill positions) plus ~20 x 1 mL QC vials. DP bag size drives LN2 sizing (open). RTx docx asks whether visual inspection can run concurrently.

RTx change and data status

MADE view

Optimization ideas

Operational impact

facility / equipment / staffing / QC

facility

equipment

staffing

QC

Data needed

Decision

Owner

Date

9 Visual inspection, cassette, controlled-rate freeze and VPLN2

C (ISO 7/8)

Day 4 into 5

Clock

RTx stated: 0.58 h*MADE working estimate: 1.75 h (rehearsal to confirm)*

Steps (RTx docx)

D4.14 · Post-fill visual inspection, label, over-wrap, vacuum seal, metal cassette, temporary fridge · 0.25 h

MADE 0.75 h (working estimate)

D4.15 · MSC close and transfer to controlled-rate freezer; CRF run; transfer to VPLN2 · 0.33 h

MADE 1 h (working estimate)

Starting material and inputs

· Labels, over-wrap, metal cassette

In-process test / sample

· Visual inspection USP <790>

Key equipment and consumable

VI light box; Vacuum sealer; Temporary fridge; Controlled-rate freezer (TF CryoMed LN2; Planer historically); VPLN2 storage

The change vs CP-1

C8 Cryopreservation: Planer CRF, current cassette to TF CryoMed LN2 CRF (or other) + metal cassette

RTx rationale: CDMO availability and broader access for pivotal. RTx strategy: Temperature mapping / equipment suitability not yet initiated; potential CDMO work package. RTx asks: How would you approach CRF mapping and validation?

Note

CRF mapping and validation not yet initiated by RTx; cassette change pending. Three-way matrix to bridge: ViaFreeze (RTx lab), LN2 CRF (EMERALD) and CryoMed (MADE).

RTx change and data status

MADE view

Optimization ideas

Operational impact

facility / equipment / staffing / QC

*facility**equipment**staffing**QC*

Data needed

Decision

Owner

Date